

PROBE MINING WITH OLIGOARRAY

SEARCHING THE DROSOPHILA HOX CLUSTERS FOR 50MER PROBE SEQUENCES

By Brian Beliveau

Last updated 8/27/12

Contact: oligopaints@genetics.med.harvard.edu

INTRODUCTION

This guide gives an overview of how to discover FISH probe sequences using the program OligoArray and assumes the user has already installed python and OligoArray along with its dependency OligoArrayAux. OligoArray and its installation instructions can be found [here](#), and OligoArrayAux and its installation instructions can be found [here](#). Instructions for downloading and installing python can be found in our [Oligopaints Scripts Manual](#). Note that this example is tailored to the LSF system employed by the Harvard RITG Orchestra UNIX cluster; users will need to modify the paths specified in this example (e.g. /files/Genetics/Wu Lab/Oligopaints/hox_example/) to match the desired working directories on their home machines or research clusters.

GETTING STARTED

The Drosophila Hox genes are broken into two ~300 kb clusters on chr3R – the antennapedia cluster which spans from *labial* to *antennapedia* and the bithorax cluster which spans from *ultrabithorax* to *abdominal B*. In dm3 coordinates, the antennapedia cluster spans from 2,487,149 - 2,824,950 and the bithorax cluster spans from 12,482,345 - 12,797,958. We can expand our target region to include flanking regions as well, such that we will mine a 1 Mb region for each target. In order to run OligoArray, we need to prepare two types of files: input files containing the sequence we want to mine for probes and a BLAST database which will be used to assess the specificity of any candidate probe sequences.


PREPARING INPUT FILES

First, we will download fasta format files for each region we want to mine from the [UCSC genome browser](#). For each cluster, we will prepare a .bed files specifying the sequences we want to download. These files will look like this:

```
track name="Antennapedia_region"  
chr3R 2156050    3156050
```

track name="biothorax_region"
chr3R 12140152 13140152

We can upload these files as Custom Tracks to UCSC:




genome: D. melanogaster assembly: Apr. 2006 (BDGP R5/dm3) [dm3]

Name	Description	Type	Doc	Items	Pos	delete
biothorax_region	biothorax_region	bed		1	chr3R:	<input type="checkbox"/>
Antennapedia_region	Antennapedia_region	bed		1	chr3R:	<input type="checkbox"/>

add custom tracks
go to genome browser
go to table browser

We can then go to the table browser and use these .bed files to download sequence files:



Use this program to retrieve the data associated with a track in text format, to calculate intersections between this form, the [User's Guide](#) for general information and sample queries, and the OpenHelix Table Browser examine the biological function of your set through annotation enrichments, send the data to [GREAT](#). Ref [Sequence and Annotation Downloads](#) page.

clade: Insect genome: D. melanogaster assembly: Apr. 2006 (BDGP R5/dm3)

group: Custom Tracks track: Antennapedia_region manage custom tracks track hubs

table: ct_AntennapediaRegion_5238 describe table schema

region: genome position chr3R:12757434-12802461 lookup define regions

identifiers (names/accessions): paste list upload list

filter: create

intersection: create


correlation: create

output format: sequence Send output to Galaxy GREAT


output file: antennapedia.fas (leave blank to keep output in browser)

file type returned: plain text gzip compressed

get output summary/statistics



Which yields a .fas file with a fasta header and our sequence of interest:



In our bioinformatics pipeline, we break our input files into 1 kb chunks using the script `input_blocks.py`:



This generates files ready to input into OligoArray.

GENERATING A CUSTOM BLAST DATABASE

Now that our input files are ready, we need to make a custom BLAST database. This step only needs to be done once for each genome to be searched. The starting point for this task is a .fas file containing the entire genome of the organism in question, arranged in the same 1 kb input file format. These can be created by downloading whole-chromosome fasta files either from UCSC using a .bed or directly from a repository such as GenBank. These files can then be concatenated using the `ConcatenateFiles.class` script available from the [OligoArray website](#). Once a master file for the whole genome is assembled, the `input_blocks.py` script can then be run. Master files for each genome already run by the Wu lab are available on the [Oligopaints website](#). In this example, a master file for the *Drosophila* genome, “dm3.fas,” already exists in our working directory. We make a new folder called “BlastDb,” in which we create a symbolic link back to our dm3.fas master file:

```

bjb11@balcony:/files/Genetics/Wu_Lab/OligoPaints/hox_example$ mkdir BlastDb
bjb11@balcony:/files/Genetics/Wu_Lab/OligoPaints/hox_example$ ls
antennapedia.fas antennapedia.txt bin biothorax.txt bithorax.fas BlastDb dm3.fas in_antennapedia.fas in_bithorax.fas
bjb11@balcony:/files/Genetics/Wu_Lab/OligoPaints/hox_example$ cd BlastDb/
bjb11@balcony:/files/Genetics/Wu_Lab/OligoPaints/hox_example/BlastDb$ ln -s ../dm3.fas
bjb11@balcony:/files/Genetics/Wu_Lab/OligoPaints/hox_example/BlastDb$ ls
dm3.fas

```

We can now make our BLAST database using the `formatdb` command with the following parameters: “-o T -p F.” Note that `formatdb` is dependent on an appropriate version of BLAST being installed. Also see the [OligoArray website](#).

```

bjb11@balcony:/files/Genetics/Wu_Lab/OligoPaints/hox_example/BlastDb$ formatdb -i dm3.fas -o T -p F
bjb11@balcony:/files/Genetics/Wu_Lab/OligoPaints/hox_example/BlastDb$ ls
dm3.fas dm3.fas.nhr dm3.fas.nin dm3.fas.nsd dm3.fas.nsi dm3.fas.nsq formatdb.log

```

RUNNING OLIGOARRAY

With our input files and BLAST database in hand, we are ready to run OligoArray. See the [OligoArray website](#) for details on all the available options. In this example, OligoArray will be run on Harvard Orchestra UNIX cluster. For each region, we will create a tab-delimited .txt file specifying the parameters to run OligoArray with. The contents of an example file, “batch_antennapedia.txt,” are show below:

```

java -Xmx8192m -jar /home/bjb11/OligoArray2_1/OligoArray2.jar -i
in_atennapedia.fas -d BlastDb/dm3.fas -o antennapedia_oligo.txt -r
antennapedia_rejected.fas -R antennapedia.log -n 30 -l 50 -L 50 -D 1000 -t 85 -T
99 -s 70 -x 70 -p 35 -P 80 -m "GGGG;CCCC;TTTTT;AAAAA" -g 52

```

This text specifies to search each region for 50mer probe sequences with an estimated T_M of between 85 and 99°C and a GC content of 35 - 80% that do not contain the sequences “GGGG,” “CCCC,” “TTTTT,” or “AAAAA.” A minimum spacing of 2 bp between probes is enforced. OligoArray will output three files: `antennapedia_oligo.txt`, which will detail every probe sequence discovered; `antennapedia_rejected.fas`, the input sequence blocks in which 0 probes were found in a format that is ready to be re-inputted into OligoArray; `antennapedia.log`, which details the line-by-line activity of OligoArray.

At this point, the files can be submitted to the cluster or run on the user’s home machine. Note that is important to make sure OligoArray2.jar and our input files have the appropriate permissions. Permissions can be checked by typing “`ls -l`,” and changed with the `chmod` command.

To run OligoArray locally, the user can simply type:

```
$ bash batch_antennapedia.txt
```

With our LSF system, the command on the cluster would be:

```
bjb11@balcony:/files/Genetics/Wu Lab/OligoPaints/hox_example$ bsub -q shared_id bash batch_antennapedia.txt
Job <14740> is submitted to queue <shared_id>.
bjb11@balcony:/files/Genetics/Wu Lab/OligoPaints/hox_example$ bsub -q shared_id bash batch_bithorax.txt
Job <14741> is submitted to queue <shared_id>.
bjb11@balcony:/files/Genetics/Wu Lab/OligoPaints/hox_example$ bjobs
JOBID USER STAT QUEUE FROM_HOST EXEC_HOST JOB_NAME SUBMIT_TIME
14740 bjb11 RUN shared_id balcony.orc clarinet001 *pedia.txt Aug 27 13:48
14741 bjb11 RUN shared_id balcony.orc clarinet001 *horax.txt Aug 27 13:48
```


Be advised that OligoArray can be very computationally intensive and can generate large (gigabyte-scale) quantities of data.

PROCESSING THE RESULTS

Once OligoArray has finished running, there were be several additional files in the working directory:

```
bjb11@balcony:/files/Genetics/Wu Lab/OligoPaints/hox_example$ ls -sh
total 286M
1.3M antennapedia.fas      2.0K antennapedia_rejected.fas  26K batch_bithorax.txt  1.4M bithorax.fas      2.0K bithorax_rejected.fas  1.4M in_antennapedia.fas
54M antennapedia.log      26K antennapedia.txt           34K bin                56M bithorax.log      34K BlastDb             1.3M in_bithorax.fas
2.4M antennapedia_oligo.txt 26K batch_antennapedia.txt     26K bithorax.txt       2.1M bithorax_oligo.txt 167M chr3.fas
```

Fortunately, the antennapedia_rejected.fas file is empty, indicating that OligoArray was able to find at least 1 probe in every 1 kb chunk inputted:



OligoArray also produced a very dense log file, a portion of which is shown below:

```
180   ...
180   Folding CCTTGGGACTCATTGCCCTAAAGACCTGGAAGGCCATGACCTTGGGCCTG... DONE
181   range=chr3R:0-999  785 rejected due to a secondary structure (deltaG = -1.459 @ 70.0 degrees)
182   Updating thermodata... DONE
183   Folding GCCTTGGGACTCATTGCCCTAAAGACCTGGAAGGCCATGACCTTGGGCCT... DONE
184   range=chr3R:0-999  784 rejected due to a secondary structure (deltaG = -1.076 @ 70.0 degrees)
185   Updating thermodata... DONE
186   Folding GGCCTTGGGACTCATTGCCCTAAAGACCTGGAAGGCCATGACCTTGGGCC... DONE
187   Testing specificity... DONE
188   range=chr3R:0-999  783 Oligo selected
189   Updating thermodata... DONE
190   Folding TTCTGCCCTTCTTCTGGGCCTGAAACTGAAGACCACTGTGCTGGTGCCA... DONE
191   Testing specificity... DONE
192   range=chr3R:0-999  731 Oligo selected
193   range=chr3R:0-999  679 rejected due to prohibited sequences
```

And finally, OligoArray outputs a file detailing the selected oligos and their properties:

