

Eric Joyce  
8/29/13

### Sameday Oligopaint production using Lambda Exonuclease

#### 1) 10 ml PCR

1 ml 10X Kapa Taq buffer  
200  $\mu$ l 10mM dNTP  
50  $\mu$ l 200uM 5'fluor-labeled Forward primer (~20-mer)  
50  $\mu$ l 200uM 5'Phosphorylated Reverse primer (~20-mer)  
10  $\mu$ l 1ng/ $\mu$ l template [Note: If using Secondaries, this will have to be a pre-touched up library]  
100  $\mu$ l Kapa Taq  
8590  $\mu$ l ddH<sub>2</sub>O  
Aliquot into strip tubes or plate  
PCR Program: 95C-5min, (95C-30s, 60C-30s, 72C-15s)X40, 72C-5min

#### 2) PCR purification by Zymo DNA clean & Concentrator-500 columns (D4031 & D4032)

Pool PCR and add 5X (50ml) DNA Binding Buffer  
Place 2 columns in 2 50ml conical tubes  
Transfer 30mls of sample to columns (15ml each) and spin 5min 4,000rpm  
Discard flow-through and add final 30ml of sample to the columns (15ml each) and spin again  
Discard flow-through and add 10ml DNA wash buffer, and spin 5 min  
Transfer columns to new tubes and Elute each with 2800  $\mu$ l ddH<sub>2</sub>O  
Should recover ~2200  $\mu$ l per column, combine for 4400  $\mu$ l (fill to this vol if necessary)

#### 3) Lambda Exonuclease Reaction to make ssDNA probe

4400  $\mu$ l DNA Eluant  
500  $\mu$ l 10X Lambda buffer  
100  $\mu$ l Lambda exonuclease 5000U/ml (if scaling up/down, keep enzyme to 5U/100  $\mu$ l of PCR)  
Aliquot into strip tubes  
Set thermo-cycler to 37C-30 min, 75C-10 min

#### 4) ppt ssDNA probe

Combine Exo reactions ~5mls  
Add 1/50 (100 $\mu$ l) glycogen and mix  
Add 1/10 (500 $\mu$ l) 4M NH<sub>4</sub>OAc and mix  
Add 2.5X (12.5mls) cold 100% EtOH and mix  
Aliquot into 2ml tubes (~9 tubes) and place at -80C for 30min  
Spin max speed at 4C for 1 hour  
Resuspend each pellet in 10 $\mu$ l ddH<sub>2</sub>O for 15min at 42C and combine  
Spec based on fluorescence for concentration  
Expect yield of ~3700pmol/10ml prep